

POWERASSIST



The PowerAssist features an easy-to-use squeeze trigger mechanism along with a speed dial, enabling users to customise speed settings.



POWER UP

The Pride Power Assist is designed to aid caregivers in pushing a manual wheelchair. The PowerAssist attaches directly on to manual wheelchairs easily.

 **FIRST
CHOICE
MOBILITY**

**027 4 CHOICE
027 424 6423**

POWERASSIST

SPECIFICATIONS

Maximum Speed¹: Up to 5 km/h (3.1 mph)

Range per Charge¹: Up to 10 km/h (6.2 mph)

Maximum Gradient: 6°

Weight Capacity: 120 kg (264 lbs)

Motor: 12V 100W 4100RPM

Total Weight: 14 kg (30.8 lbs)

Battery: 12V 17ah Sealed Lead Acid

Battery Weight: 6.05 kg (13.3 lbs)

Charger: 3A 12V output


Fits Wheelchair Widths: 16" - 20"

Dimensions ² :	Height:	285mm (11.2")
	Length:	280mm (11")
	Width:	220mm (8.6")

¹ Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure.

² Due to manufacturing tolerances and continued product improvement, this specification is subject to a variance of +/- 3%.

The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development, they reserve the right to alter specifications without prior notice.



The Power Assist Battery Pack



The Power Assist Unit

FEATURES

- The Pride Power Assist is designed to aid caregivers in pushing a manual wheelchair in otherwise strenuous situations such as up steep hills, inclines and loose gravel.
- The Power Assist is designed to work in conjunction with most wheelchairs with seat widths of 16"-20".
- Weight capacity of 120 kg/264 lbs (wheelchair and user combined).
- The Power Assist has a maximum speed of up to 5 km/h (3.1 mph).

 **FIRST
CHOICE
MOBILITY**
027 4 CHOICE
027 424 6423