AML Convertible Rollator Transfer Chair





Code: AMLTRROL

Description

This lightweight collapsible aluminium rollator has four 8.0" castors. The front two castors swivel allowing for the rollator to be easily moved in any direction. These larger sized castors allow it to be used both indoors and outdoors.

The cable loop brakes are easy to use and when squeezed with both hands allow instant braking. They also have a locking feature which allows the brakes to be "locked on", necessary when using the wheeled walker's built-in seat.

This rollator features a transfer chair with a height adjustable footrest and padded backrest.

Features

High strength aluminium tubing Ergonomic Hand Grips Height adjustable Footrest Transfer chair Padded seat and backrest Under seat storage Locking looped hand brakes Folds up for ease of transport Some assembly required



09 236 4233 027 424 6423 / 027 4 CHOICE firstcmobility@gmail.com www.firstchoicemobility.co.nz

Specifications

Recommended Maximum User Weight 136 kg
Seat Height
Seat Size (W x D) 346.5 x 350 mm
Handle Height Adjustment 825 - 975 mm
Overall Width 620 mm
Boxed Dimensions (W x D x H) 215 x 630x 885 mm
Unit Weight 8.7 kg
Colours Available Blue