

Pride® QUEST

FOLDING SCOOTER

THE **ULTIMATE**
COMPACT, **FOLDABLE**
TRAVEL **SCOOTER**

Ready for transport and easily folded, the 4-wheel **Quest Folding Scooter** is perfect for the active lifestyle. Offering great performance indoors and out at a maximum speed of 6 km/h (4 mph) and a range per charge of up to 20 km (12.5 miles), the Quest Folding Scooter provides a comfortable and convenient ride. It is super portable and is a great space saver that adjusts to fit perfectly inside any closet, corner, or vehicle trunk. The Quest Folding Scooter is sleek, stylish, and extremely easy to operate, making it the perfect travel mobility scooter choice.



Scooter Folded



Easy to Transport



**FIRST
CHOICE
MOBILITY**
027 4 CHOICE
027 424 6423

PRIDE® **MOBILITY** SCOOTERS
LIVE YOUR BEST®

SPECIFICATIONS

QUEST FOLDING SCOOTER

4W: S19

| | |
|---|--|
| Weight capacity: | 115 kg (250 lbs.) |
| Maximum speed ¹ : | Up to 6 km/h (4 mph) |
| Ground clearance ² : | 6 cm (2.3") |
| Turning radius ² : | 82 cm (32") |
| Overall length ² : | 93 cm (37") |
| Overall width ² : | 48.5 cm (19") |
| Drive wheels ² : | 5 cm x 20 cm (2" x 8") solid |
| Front wheels ² : | 4 cm x 18 cm (1.6" x 7") solid |
| Range per charge ^{1,4} : | Up to 20 km (12.5 miles) (lithium-ion); |
| Total weight without batteries ⁶ : | 23.5 kg (51.4 lbs.) |
| Heaviest piece when disassembled: | 23.5 kg (51.4 lbs.) (unit does not disassemble - unit weight without battery pack) |
| Standard seating: | Type: Foldable, compact, vinyl-covered molded plastic Material: Gray Vinyl Dimensions: 41 cm x 33 cm (16" x 13") |
| Battery requirements ⁵ : | (1) 24-volt (lithium-ion) Size: (1) 11.5 AH (lithium-ion) Weight ³ : 11.5 AH lithium-ion- 2.3 kg (5 lbs.) each (1 required) |
| Battery charger: | 2.5-amp, off-board |

Colour:



White

- ¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
² Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+/- 3%).
³ AGM or gel cell type required.
⁴ Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance test conducted at maximum weight capacity.
⁶ Battery weight may vary based on manufacturer.

All specifications provided by Heartway® Medical Products.
 Product distributed by Pride Mobility®.

NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tire pressures, ambient temperature, and the gradient and surface of the road or pavement. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

FEATURES

- Unique folding design
- Folds in just three easy steps
- Adjustable tiller
- Up to 20 km (12.5 miles) range per charge
- Easy to transport at home, hotels, or airports
- Can be stowed upright or laid down
- Articulated front end adds greater stability
- Front and rear LED lights

TO FOLD:



Fold Seat



Pull Up Lever to Fold



Fold and Lock Tiller

TO UNFOLD:



Pull Lever



Unfold Scooter



Unfold Seat; Lock in Place



Note: Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your Go-Go® Folding Scooter and its lithium-ion batteries are permitted.

FIRST CHOICE MOBILITY

BRUCE & SHARON ENGLAND

Office: 09 236 4233

Mobile: 027 424 6423 / 027 4 CHOICE

Email: firstcmobility@gmail.com

Website: www.firstchoicemobility.co.nz